thermoscientific

PRODUCT SPECIFICATIONS

Thermo Scientific 49iQPS

UV Photometric Primary Standard—calibration of ozone analyzers and transfer standards



The Thermo Scientific[™] 49iQPS Ozone (O₃)Primary Standard utilizes UV Photometric technology to produce ozone at a rate of up to 5000 ppb.

The Thermo Scientific 49iQPS
Analyzer is a dual cell photometer,
the concept adopted by NIST for
the national ozone standard. The
49iQPS Analyzer can operate with
ozonator flow rates of up to 6 liters a
minute. Because the instrument has
both sample and reference flowing
simultaneously, a response time
of 20 seconds can be achieved.
Temperature and pressure correction
are standard offerings.



Non-Stop Intelligence

- Predictive Diagnostics
- Proactive Communication
- Personal Device Connectivity

The Thermo Scientific iQ Series Gas Analyzer provides a smart environmental monitoring solution designed for reliability, easy operation and proactive maintenance. Get more control over your instrument's performance, costs, workflow and data availability.



The iQ companion app for the iQ Series Gas Analyzer delivers the ultimate in ease of use and smart engineering. The iQ app allows for remote monitoring of iQ gas analyzers, simplified

ways of contacting us and instant access to product resources. Download the iQ app at **thermofisher.com/iQapp**

Thermo Fisher

thermoscientific

Thermo Scientific 49iQPS Ozone Primary Standard

	orar o ozono i rimar y otanicar a
Specifications	
Range	50-5000 ppb; 100-10000 25 μg/m³
Zero noise	0.25 ppb RMS (60 second averaging time)
Detection limit	1.0 ppb (60 second averaging time)
Response time	20 seconds (0-95%)
Lag time	10 seconds
Linearity	±1% full scale
Flow rate	1-3 LPM
Operating temperature	0°C - 45°C
Power requirements	100-240 VAC 50/60 275 Watts
Size and weight	16.75" (W) \times 8.72" (H) \times 24" (D), 34.53 lbs. 425.25 mm (W) \times 221.48 mm (H) \times 609 mm (D), 15.7kg
Analog I/O	4 Isolated voltage inputs 0–10 V 6 Isolated analog voltages outputs, with 4 selectable ranges 6 Isolated analog current outputs, with 2 selectable ranges
Digital I/O	16 Digital inputs (TTL)8 Solenoid driver outputs10 Digital reed relay contact outputs
Serial ports	1 RS-232/485 port; 1 RS-485 external accessory port
Other ports	3 Full Speed USB ports (one in front, two in rear) 1 Gigabit ethernet port
Communication	MODBUS, Streaming, AK, Gesytec (Bayern-Hessen)
Approvals and Certifications	CE, TUV-SUD Safety
Ozonator	
Output	.025 - 1.000 ppm @ 3-4 LPM
Response	1 minute to 98% or 5 ppb of final value, whichever is greater
Stability	±4 ppb or ±1% of reading, whichever is greater

Ordering information

49iQPS Ozone Primary Standard

Choose from the following configurations/options to customize your own 49iQPS Analyzer

1. Power Cord				
A = 100-120 VAC 50/60 Hz (NA)				
B = 220 VAC 50/60 Hz (CHN)				
C = 220 VAC 50/60 Hz (EU)				
2. Communications				
N = No I/O				
A = Serial RS232/RS485				
B = Analog and Digital				
C = Serial, Analog and Digital				
3. Zero air source:				
NN = No Zero Air Source				

1	1	l

NA = Zero air source

Your Order Code: Model 49iQPS -

To maintain optimal product performance, you need immediate access to experts worldwide, as well as priority status when your air quality equipment needs repair or replacement. We offer comprehensive, flexible support solutions for all phases of the product life cycle. Through predictable, fixed-cost pricing, our services help protect the return on investment and total cost of ownership of your Thermo Scientific products. For more information on our comprehensive service solutions visit thermofisher.com/EMservice

USA27 Forge Parkway
Franklin, MA 02038
Ph: (508) 520-0430
Fax: (508) 520-2800

orders.aqi@thermofisher.com

India
C/327, TTC Industrial Area
MIDC Pawane
New Mumbai 400 705, India
Ph: +91 22 4157 8800
india@thermofisher.com

China

+Units 702-715, 7th Floor Tower West, Yonghe Beijing, China 100007 Ph: +86 10 84193588 info.eid.china@thermofisher.com Europe

Ion Path, Road Three, Winsford, Cheshire CW73GA UK Ph: +44 1606 548700 Fax: +44 1606 548711 sales.epm.uk@thermofisher.com

Find out more at thermofisher.com/49iQPS thermofisher.com/iQSeries

